


NEW!

OrthoSta

Ligament Repair System

The versatile ligament repair system for stabilising injured joints.

- High-strength woven UHMWPE fiber tape available in 2 sizes
- Titanium interference screw fixation for knotless and crimpless system
- Features Healix™ antimicrobial technology for infection prevention



Products of the OrthoSta Ligament Repair System

PATIENT WEIGHT RECOMMENDATION < 18 Kg — OrthoSta System

Item#	Description
OSFTL 2036	<i>OrthoSta UHMWPE Woven Tape 2mm, Blunt Needle, 36cm, TI SB-AS, Healix Antimicrobial</i>
OSIS 200.08	OrthoSta Interference Screw, for 2mm Tape, 8mm, T6 Star Recess, TI-6AL-4V
OSTAP 2.0	OrthoSta Tap, for 2mm Tape, Stainless Steel
500035	Screwdriver Blade, Star Driver QCK Shaft, T6
2.5 QCK 110	2.5mm Quick Coupling Drill Bit, Tin Coated, 110mm - Material 440A

PATIENT WEIGHT RECOMMENDATION > 18 Kg — OrthoSta System

Item#	Description
OSFTL 4060	<i>OrthoSta UHMWPE Woven Tape 4mm, Blunt Needle, 60cm, TI SB-PS, Healix Antimicrobial</i>
OSIS 400.13	OrthoSta Interference Screw, for 4mm Tape, 13mm, T8 Star Recess, TI-6AL-4V
OSTAP 4.0	OrthoSta Tap, for 4mm Tape, Stainless Steel
500036	Screwdriver Blade, Star Driver QCK Shaft, T8
3.5 QCK 110	3.5mm Quick Coupling Drill Bit, 110mm - Material 440A



Learn to Use the OrthoSta Ligament Repair System

Visit education.movora.com for upcoming workshops featuring the OrthoSta system! Start learning the system today with OrthoSta Knotless Ligament Repair System Online Course



Scan Here

To View OrthoSta
Online Course

